




INSTI[®]
HCV Antibody Test

Introducing the World's First One Minute
HCV Antibody Test

- + Easy to use with results in one minute
- + Clear and accurate results
- + Enhances patient flow
- + Minimizes face to face contact time



0543






HCV Antibody Test

Detects HCV genotypes 1 - 6


IgG/IgM Procedural Control

All support materials included


Sample & pour



Read results immediately!



Non-Reactive
1 Control dot



Reactive
1 Control dot
1 Test dot

Product Information

Information Type	Product Details
Method	Flow-through
Time to results	As little as 60 seconds
Storage conditions	35.6–86 °F (2–30 °C)
Shelf life	12 months
Sample type	Fingerstick whole blood, venipuncture whole blood, plasma or serum

Accuracy

Sample Type	Sensitivity	Specificity
Fingerstick Whole Blood	100%	100%
Serum	100%	99.4%
Plasma	100%	99.6%

Ordering Information

Product

90-1062-001 - INSTI® HCV 25 Antibody Test Kits

90-1062-002 - INSTI® HCV 50 Antibody Test Kits

90-1062-003 - INSTI® HCV 300 Antibody Test Kits

90-1105 - INSTI® HCV Bulk Components 50-Pack

Manufactured by:



bioLytical Laboratories Inc.

406-13251 Delf Place

Richmond, BC Canada V6V 2A2

1-866-674-6784 | info@biolytical.com

insti.com



65-5001(D) 04/2023