




Multiplex
Syphilis
Ab Test

HIV-1/2



Syphilis

Control

HIV

Introducing The World's Fastest
Dual Point-of-Care Test

- + Superior speed
- + Greater workflow flexibility
- + 1 sample, 1 minute, 2 results



Is 60 seconds quick enough?

insti.com • biolytical.com

INSTI® Multiplex HIV-1/2 Syphilis Ab Test

Unique “flow-through” design


IgG Capture Procedural Control

No timers required


Measures antibodies for *T. pallidum*, in addition to HIV antibodies

Knowing your result can be as easy as 1, 2, 3...

Sample & pour



Read results immediately!



Product Information

Information Type	Product Details
Method	Flow-through
Time to results	As little as 60 seconds
Storage conditions	15–30 °C (59–86 °F)
Shelf life	12 months
Sample type	Fingerstick, venous whole blood, plasma or serum

Accuracy¹

	Sensitivity	Specificity
HIV-1/HIV-2	100%	99.5%
Syphilis ²	99.5%	99.9%



Confidence. The built-in control technology confirms correct sample addition, test procedure and proper test performance.

Ordering Information

Product

INSTI® Multiplex HIV-1/2 Syphilis Ab Test

INSTI® HIV-1/HIV-2 Positive and Negative Controls (3 x 1.0mL vials)

INSTI® *T. pallidum* Control (1 x 1.0mL vial)

For ordering information, please contact:

+1 604 204 6784
info@biolytical.com

References

¹INSTI Multiplex HIV-1/2 Syphilis Ab Test Package Insert (CE), 51-1143 (K).

²In cases of active syphilis when RPR>1:4